



DETECTO's model 450 baby scale shown here with 130 lb x 1 oz capacity and built-in measuring tape.



450 SERIES

MECHANICAL BABY SCALES

- Time-tested, reliable weighbeam accuracy
- No power required – ideal for remote clinics
- Baked-on, powder painted steel base and column
- Comfortable, polystyrene weighing tray
- Built-in measuring tape:
0 to 21 in x 0.25 in /
0 to 53 cm x 0.1 cm



BUILT-IN MEASURING TAPE

The easy-to-read measuring tape allows users to measure both weight and height simultaneously (0 - 21 in/53 cm).



TIME-TESTED WEIGHBEAM

DETECTO's die-cast mechanical weighbeam is world-renowned for its high accuracy and long-term dependability over decades of use.



READABLE FROM BOTH SIDES

The weighbeam is imprinted on both sides, so the nurse can read the baby's weight from either side, front or back.

MECHANICAL BABY SCALES

DURABLE CONSTRUCTION / CLINICAL QUALITY

Classic design and high-quality precision mechanics

DETECTO's mechanical baby scales provide the ultimate in safety for pediatric patients while also utilizing the time-tested, reliable DETECTO weigh beam for highly-accurate measurement. They are available in capacities up to 130 pounds and 65 kilograms with your choice of pounds or kilograms readings. There is no electricity required, making them ideal for remote clinics and pediatrician offices.

The easy-to-clean polystyrene tray cradles babies securely on these dependable pediatric scales. The heavy-duty steel base provides patient assurance while weighing. A built-in measuring tape allows the healthcare professional to measure both weight and height simultaneously.

The die-cast weigh beam is positioned high above the scale base to provide for unobstructed viewing of the baby's weight measurement at an easy-to-read height. The baked-on, white powder paint offers a long-lasting dependable finish. For fast mechanical measurement, DETECTO is the most trusted name brand in clinical weighing.

Model	450	451
Capacity	130 lb x 1 oz	65 kg x 20 g
Measuring Tape	0 to 21 in x 0.25 in / 0 to 53 cm x 0.1 cm	
Dimensions	Overall	22 in W x 21 in D x 27 in H / 56 cm W x 53 cm D x 69 cm H
	Tray	22 in W x 14.75 in D x 4.5 in H / 56 cm W x 37 cm D x 11 cm H
	Platform	10.75 in W x 14.75 in D / 27 cm W x 37 cm D
Construction	Base and Column	White baked-on, powder-painted steel
	Weighing Tray	White polystyrene
Net Weight	29 lb / 13 kg	
Shipping Weight	33 lb / 15 kg	
UPC Code	809161142402	809161142501

DETECTO reserves the right to improve, enhance, or modify features and specifications without prior notice.



ROOMY WEIGHING TRAY

The comfortable polystyrene tray cradles babies securely while being weighed and measured.



 **DETECTO**
www.Detecto.com

102 E. Daugherty, Webb City, MO 64870 USA
Ph: 417-673-4631 or 1-800-641-2008 Fax: 417-673-2153
www.Detecto.com

SOLD BY: