DETECTO 6100



PORTABLE STAND-ON SCALES



PORTABLE PATIENT SCALE WITH CARRYING HANDLE AND REMOTE INDICATOR

FEATURES

- 800 lb x 0.2 lb / 360 kg x 0.1 kg capacity
- 18 in W x 14 in D x 2 in H / 46 cm W x 36 cm D x 6 cm H platform
- · Integral carrying handle for portability
- Body Mass Index calculation
- -C connected models include Wi-Fi and Bluetooth for wireless connectivity
- · -AC models include AC adapter
- · Four models available to choose from
- · Bubble level helps you find level weighing foundation for optimal accuracy
- · Non-skid black mat for patient safety
- Lock/Release key can hold weights and BMI calculation
- Indicator may be situated on a nearby tabletop, counter, or floor
- · Large, easy-to-read 0.75-in/19-mm high LCD weight digits





102 E. Daugherty, Webb City, MO 64870 USA Ph: 417-673-4631 or 1-800-641-2008 Fax: 417-673-2153 **SOLD BY:**