

DETECTO 6100

PORTABLE STAND-ON SCALES

DETECTO®



DETECTO 6100

PORTABLE PATIENT SCALE WITH CARRYING HANDLE AND REMOTE INDICATOR

FEATURES

- 800 lb x 0.2 lb / 360 kg x 0.1 kg capacity
- 18 in W x 14 in D x 2 in H / 46 cm W x 36 cm D x 6 cm H platform
- Integral carrying handle for portability
- Body Mass Index calculation
- -C connected models include Wi-Fi and Bluetooth for wireless connectivity
- -AC models include AC adapter
- Four models available to choose from
- Bubble level helps you find level weighing foundation for optimal accuracy
- Non-skid black mat for patient safety
- Lock/Release key can hold weights and BMI calculation
- Indicator may be situated on a nearby tabletop, counter, or floor
- Large, easy-to-read 0.75-in/19-mm high LCD weight digits



Powered by 6 AA batteries
making the scale ultra-portable

Versatile indicator mounting
bracket included for portable
countertop or permanent wall
mounting use

DETECTO®

102 East Daugherty St., Webb City, MO 64870 USA
Toll-Free: (800) 641-2008 E-mail: detecto@cardet.com
www.Detecto.com

SOLD BY: