

# 6856

## Bariatric Waist-High Portable Medical Scale

Designed for obese and unsteady patients, this 1,000-lb capacity scale features a balanced, lightweight design. It's easily tipped back for transportation on its hard, large diameter, rubber wheels yet narrow enough to fit through standard doors. The battery-powered, LCD display features easy-to-use Body Mass Index calculation while weighing. The low-profile platform, convenient handrail, sturdy steel construction, and skid-proof mat add to patient assurance and stability.



### SPECIFICATIONS

<b>Model:</b>	6856
<b>Capacity:</b>	1,000 lb x 0.2 lb / 450 kg x 0.1 kg
<b>Platform Size:</b>	24"W x 24"D / 61 cm x 61 cm
<b>Display:</b>	<ul style="list-style-type: none"><li>• 0.75" high digits, high-contrast LCD</li><li>• Body Mass Index calculation</li><li>• LOCK button holds weight after weighing</li></ul>
<b>Power:</b>	6 "AA" size batteries (not included) or optional AC adapter.
<b>Shipping Weight:</b>	85 lb / 39 kg

Options include height rod and AC adapter.



*Large LCD and Body Mass Index calculation*