

OIML EC Class III model also available (6868KGEU)

E			588 FRICE LEE	
	Capacity	1,000 lb x 0.2 lb / 450 kg x 0.1 kg	300 kg x 0.1 kg or 350 kg x 0.2 kg (OIML EC Class III)	
SIONS	Overall Dimensions			
DIMENSIONS	Platform Dimensions			
	Construction	Aluminum and stainless steel		
	Power	6 AA batteries (not included) or optional AC adapter	6 C batteries (not included) or optional AC adapter	
OR	Display	Extra-large, 0.75 in/ 19 mm-high LCD digits		
INDICATOR	Keypad	Membrane type with 5 keys and 4 directional arrows		
IND	Buttons	Zero, Units, Lock/Release, On/Off, BMI/Enter, Print, Arrow keys		Rubber mat provides
	Connectivity	RS232 serial port		traction and stability
	-	80 lb / 36 kg		
VEIGHT	Net Weight	80 lb / 36 kg		(4) load cells



The 6868 bariatric scale uses 6 AA batteries (not included) or an optional AC adapter (model 6800-1045). DETECTO's model 6868 clinical scale offers a variety of usages due to its unique product design and high capacity: standing, seated, and bariatric weighing. The wraparound stainless steel handrails with padded grips provide comfortable support for obese and unsteady patients. The 6868 is completely portable and can be transported easily due to its battery power and integral wheels (6 AA batteries not included).

- Body Mass Index calculation
- Lock/Release key for viewing weight after the patient steps off the scale
- Automatic Shutoff and Sleep modes conserve battery life
- Sturdy, long-lasting quality construction
- Many other bariatric scales also available from DETECTO
- Fits through doorways up to 32 inch (81 cm) wide



203 E. Daugherty, Webb City, MO 64870 USA Ph: 417-673-4631 or 1-800-641-2008 Fax: 417-673-2153

Mobile: m.detecto.com

DETECTO reserves the right to improve, enhance, or modify features and specifications without prior notice.

**SOLD BY:**