

DETECTO 750 Clinical Weight Indicator

SERIAL CONNECTIVITY STANDARD

Serial ports on all model 750 indicators used with DETECTO medical scales.



Serial I/O Connector

RS232 serial port allows for electronic medical recordkeeping



DETECTO's model 750 digital weight indicators provide RS232 serial connectivity as a standard feature. This allows any DETECTO medical scale utilizing the 750 indicator to have EMR/EHR-ready serial output for interfacing with a computer, monitor, or printer for efficient data transfer in electronic medical recordkeeping.



**Scale
Connectivity for
Patient Records**

- Body Mass Index calculation
- Extra-large, 0.75 in/19 mm high LCD digits
- 6 AA batteries (not included) or optional AC adapter
- Lock/Release key for viewing weight after the patient steps off the scale
- Units conversion for lb or kg
- Durable steel construction for years of reliable service

