


84435 DIGITAL BABY SCALES

DIGITAL SPEED AND ACCURACY





HEAVY-DUTY SOLID, STEEL SCALE BASE

for patient assurance while weighing.

- Pushbutton tare to remove blanket weights
- Built-in measuring tape: 0 to 21 in x 0.25 in / 0 to 53 cm x 0.1 cm
- Easy-to-clean plastic trays for baby comfort
- Easy-to-read highcontrast LCD
- Automatic weight hold feature for active babies, locks in weight
- Heavy steel base for patient stability
- Durable, long-lasting steel construction
- AC or battery power

DETECTO DIGITAL BABY SCALES

DURABLE CONSTRUCTION / CLINICAL QUALITY



Model	8435
Capacity	130 lb x 0.1 lb / 59 kg x 0.05 kg
Overall Dimensions	37.5 in H x 22 in W x 21 in D / 95 cm H x 56 cm W x 53 cm D
Measuring Tape	0 to 21 in x 0.25 in / 0 to 53 cm x 0.1 cm
Tray Size	22 in L x 14.7 in W x 4.5 in D (56 cm L x 37 cm W x 11 cm D)
Display	Six-digit, seven segment, 0.7 inch (17.8 mm) high LCD
Power	Six "AA" Alkaline, Ni-Cad or NiMH batteries (not included) OR an optional 100 to 240V AC 50/60Hz 12V DC 1A wall plug-in UL/CSA listed AC power adapter (Cardinal part number 6800-1045).
Keyboard	Membrane type with 5 keys and 4 directional arrows
Connectivity	RS232 serial port
Construction	Baked-on, powder-painted steel
Shipping Weight	33 lb (15 kg)





height simultaneously.

Detecto's model 8435 features an easy-to-use keypad and extra-large display perfect for patient weighing and viewing by the attending nurse, physician, or home healthcare professional.



CONVENIENT DISPLAY POSITIONING

The weight indicator is conveniently situated over the top of the child, so it's easy to read the weight while the baby is on the scale.



QUALITY CONSTRUCTION

The 8435 has a lock-in weight feature which compensates for typical baby movement. The crisp LCD features 0.75 inch high digits for easy readouts.



BUILT-IN MEASURING TAPE

The easy-to-read measuring tape allows users to measure both weight and height simultaneously (0 - 21 in/53 cm).

DETECTO reserves the right to improve, enhance, or modify features and specifications without prior notice.

DETECTO www.Detecto.com

102 E. Daugherty, Webb City, MO 64870 USA Ph: 417-673-4631 or 1-800-641-2008 Fax: 417-673-2153 www.Detecto.com

SOLD BY: