





Large, easy-to-read display with weight recall



The tray collapses for easier transport and to allow for convenient standing toddler weighing.



- Storage of previous weight memory recall
- Calculation of change in weight
- Low battery indicator
- Optional carrying case available (8440-CASE)





Date and time displays automatically in sleep mode or by pressing the Time/Weight button during weighing

8440 BABY SCALE SPECIFICATIONS

CAPACITY: 44 lb x 1/2 oz / 20 kg x 10 g

DISPLAY: 1 in / 25 mm high LCD

TRAY DIMENSIONS: 20.6 in W x 10.5 in D x 4 in H /

523.9 mm W x 266.7 mm D x 101.6 mm H

SCALE DIMENSIONS: 9.8 in W x 11.3 in D x 7.5 in H /

247.7 mm W x 285.8 mm D x 190.5 mm H

SHIPPING WEIGHT: 8 lb / 3.6 kg

POWER: (1) 9V Alkaline battery (not included) and (1) 3V

CR2032 lithium back-up battery (included) for memory

CLOCK FUNCTION: Time and date

MEMORY FUNCTION: Memory for comparison of current and stored weight

Detecto reserves the right to improve, enhance, or modify features and specifications without prior notice.

A Division of Cardinal Scale Manufacturing Co. 203 E. Daugherty, Webb City, MO 64870 USA

Ph: 417-673-4631 or 1-800-641-2008 Fax: 417-673-5001 www.Detecto.com

Sold By: