



BARIATRIC SCALES

featuring
DIGITAL
Height Rods



**DIGITAL SPEED
AND ACCURACY**

Heavy-duty, yet lightweight, the 6857DHR bariatric scale's balanced design and integral wheels make it easy to move. The 6857DHR is engineered specifically for weighing obese and unsteady patients and features a roomy, spacious 2-ft-square platform.

DETECTO
Model 6857DHR

DETECTO DIGITAL BARIATRIC SCALES

-- WITH DIGITAL HEIGHT RODS --



DIGITAL WEIGHT • HEIGHT • BMI

The 6857DHR and 6854DHR bariatric scales are designed for optimum performance and display weight, height, and BMI simultaneously (shown below) by

VERSATILE CONNECTIVITY

Detecto's state-of-the-art MedVue® digital weight indicator includes DB9 serial, USB, and optional wired or Wi-Fi Ethernet output for data transfer to a printer or PC-based EMR (electronic medical record) software. With the time/date feature, the scale can send the weight, time, date, BMI (Body Mass Index), and height directly to a printer or a PC for record keeping. The 350-degree swivel mounting bracket allows it to be positioned wherever it's easiest for viewing.

simply raising the digital height rod while a patient is being weighed on the scale. There are no keys to push, no operational sequences to memorize, it's fully automatic. These bariatric scales are EMR-ready with USB, serial, and optional Wi-Fi. The ability to print or save weight data with these scales is a real time-saving device for clinics, doctor's offices, long-term care facilities, and other medical institutions.



- ▶ Digital measurement speed and accuracy
- ▶ Integral wheels for quick and easy transport
- ▶ Stainless steel handrails for assurance
- ▶ Height rod measures up to 6'7" / 200 cm

DETECTO DIGITAL BARIATRIC SCALE SPECIFICATIONS

	6857DHR	6854DHR
CAPACITIES:	1,000 lb x 0.2 lb / 450 kg x 0.1 kg	600 lb x 0.2 lb / 270 kg x 0.1 kg
PLATFORM SIZE:	24" x 24" / 61 cm x 61 cm	18" x 14" / 46 cm x 36 cm
SHIPPING WEIGHT:	85 lb / 39 kg	74 lb / 34 kg
DIGITAL HEIGHT ROD MEASUREMENTS:	43" / 110 cm - 79" / 200 cm	same
POWER:	6 AA size alkaline batteries (not included) or optional AC adapter	same
DISPLAY:	Liquid Crystal, 0.875"/22 mm high weight readouts	same
CONNECTIVITY:	Serial, USB, and optional Wi-Fi	same



DETECTO
Model 6857DHR



DETECTO
Model 6854DHR

Detecto Scale reserves the right to improve, enhance, or modify features and specifications without prior notice.



Sold By:



102 East Daugherty St., Webb City, MO 64870 USA
Toll-Free: (800) 641-2008 E-mail: detecto@cardet.com
www.Detecto.com