

# DIGITAL STADIOMETER model DHRWM

(ProDoc)

### Standalone Wall-Mount Digital Height Rod

Detecto's exclusive digital height technology. Receive height readings quickly, accurately, and comfortably.

Measure up to 6'7" / 200 cm

Automatic height viewing on the remote display, positioned where you need it with Detecto's wall-mount or desktop bracket.





# PRODOC

## DIGITAL STADIOMETER model DHRWM

## Standalone Wall-Mount Digital Height Rod

### Features

- Fast, accurate digital height LCD readouts
- Modern, clean design
- Display can be mounted on wall or desk
- View height in inches or centimeters
- Measures up to 6'7" / 200 cm

### **Specifications:**

DISPLAY: Five-digit, seven-segment, 1 in (25 mm) high LCD MEASURES: 43 in / 110 cm - 79 in / 200 cm HEIGHT ROD INCREMENTS: 0.5 in / 1 cm POWER: 6 "AA" Alkaline batteries (not included) CABLE LENGTH: 69 in / 175 cm (display to height rod) SHIPPING WEIGHT: 8.2 lb / 3.7 kg

Detecto Scale reserves the right to improve, enhance, or modify features and specifications without prior notice.



Sold By: