

DETECTO®

ProDoc®

DIGITAL STADIOMETER

model **DHRWM**

**Standalone Wall-Mount
Digital Height Rod**

Detecto's exclusive digital height technology. Receive height readings quickly, accurately, and comfortably.

Measure up to 6'7" / 200 cm

Automatic height viewing
on the remote display,
positioned where you need
it with Detecto's wall-mount
or desktop bracket.



ProDoc DETECTO

5'09.5"

UNITS

ON I
OFF O



ProDoc[®]

DIGITAL STADIOMETER model **DHRWM**

Standalone Wall-Mount Digital Height Rod

Features:

- Fast, accurate digital height LCD readouts
- Modern, clean design
- Display can be mounted on wall or desk
- View height in inches or centimeters
- Measures up to 6'7" / 200 cm

Specifications:

DISPLAY: Five-digit, seven-segment, 1 in (25 mm) high LCD

MEASURES: 43 in / 110 cm - 79 in / 200 cm

POWER: 6 "AA" Alkaline batteries (not included)

CABLE LENGTH: 69 in / 175 cm (*display to height rod*)

SHIPPING WEIGHT: 8.2 lb / 3.7 kg

Detecto Scale reserves the right to improve, enhance, or modify features and specifications without prior notice.

Sold By:



DETECTO[®]

102 East Daugherty St., Webb City, MO 64870 USA
Toll-Free: (800) 641-2008 E-mail: detecto@cardet.com
www.Detecto.com