

# DIALYSIS PRODUCTS

High-capacity scales with innovative features designed for dialysis patients



**DETECTO**

# WIDE VARIETY OF SCALES FOR DIALYSIS PATIENT WEIGHING

## Offering the ultimate in patient satisfaction and accuracy

With dialysis scale use on the rise, DETECTO has developed state-of-the-art stainless steel scale platforms and software specially designed for dialysis patient care. Patients that are ambulatory, using a walker or cane, and patients in wheelchairs and beds can easily access the weight platform that is installed flush with the floor surface or low profile to the ground. The in-floor design is flush-mounted to the floor, and the indicators are placed within the facility's wall to save valuable space for an integral, streamlined patient weighing operation.

### ID-7248SH-855RMP

#### Solace In-Floor Flush-Mounted Dialysis Scale

Solace series in-floor dialysis scales are flush-mounted platform scales that feature stainless steel decks and include a pit frame, trim ring, and interface cable to the included multi-color, touchscreen LCD indicator and printer. The flush surface eliminates tripping hazards and hallway accessibility issues. The Solace series is available with or without a handrail and all models include a recessed-mount indicator and ticket printer combination. Special dialysis weighing software is included with the scale. These in-floor scales with high 1,000 lb x 0.2 lb / 450 kg x 0.1 kg capacity and a variety of sizes from 3 x 3 to 6 x 4 feet.



### 7550

#### Wall-Mount Fold-Up Scale

The 7550 wall-mount wheelchair scale features a fold-up design that saves space when not in use. It can be installed in high-traffic areas for convenient, weighing, and fold up when finished weighing. The scale accommodates most wheelchairs, and yet requires very little floor space when folded up. A large, colorful, interactive touchscreen allows you to enter patient IDs that can be sent along with patient weight data to EMR systems via Wi-Fi or Bluetooth 5.0. The scale also offers wheelchair tare weight mode to remove wheelchair weights and factor patient weight (pushbutton and keypad tare).



### 6856MHR

#### Digital Handrail Scale with Height Rod

The 6856MHR bariatric scale is designed for optimum weighing performance and stability with its wraparound tubular stainless steel handrails and padded handgrips. The spacious, sturdy steel platform measures 24 in/61 cm square. Integral wheels and battery power allow easy transport. Scale capacity is 1,000 lb x 0.2 lb / 450 kg x 0.1 kg. Height rod measures from 24-84 in / 60-212 cm.



### 6855

#### Compact Floor Scale with Secure Wraparound Handrails

DETECTO's digital handrail scale is made specifically for accommodating large and unsteady patients. The scale's low-profile platform features a non-skid rubber mat surface, and the scale's wrap-around stainless steel handrails come with padded grips for patient confidence. An RS232 serial port comes standard, with Wi-Fi/Bluetooth connectivity available. The scale's capacity is 600 lb x 0.2 lb / 270 kg x 0.1 kg and platform size is 18 x 14 in.



### 6475

#### Digital Chair Scale with 400 lb Capacity

The ultimate in patient safety and comfort, DETECTO's digital chair scales are completely mobile, with built-in wheels and battery power. Lift-away armrests and a fold-up footrest make it easy for transferring dialysis patients, disabled patients, and the elderly. The chair scale comes with an RS232 serial port, with Wi-Fi/Bluetooth connectivity available.



### P-48R

#### Step-On Cans, the Rugged Waste Receptacle

Named for their rugged toughness, Step-On Cans are perfect for laundry departments and receiving areas. Equipped with a heavy-duty foot pedal and linkage and a bag securing mechanism, they come in a variety of five different sizes and capacities (4, 6, 8, 12 and 25 gallon) and two different materials: stainless steel or a baked epoxy finish (available in red or white). The cans are flame retardant, have a self-closing lid for sanitation, and feature unitized construction.



### GH3

#### Universal Wall-Mount Steel Glove Box Holder

Conserve counter space with DETECTO's universal glove box holders. The open design allows you to see what size or type of glove is being dispensed. The holders accommodate most glove boxes and are universal enough to be used for tissue holders where counter or shelf space is limited. Models capable of holding one, two, or three boxes of gloves. Lifetime warranty.



*DETECTO reserves the right to improve, enhance, or modify features and specifications without prior notice.*

## SOLD BY:



102 East Daugherty St., Webb City, MO 64870 USA  
Toll-Free: (800) 641-2008 E-mail: [detecto@cardet.com](mailto:detecto@cardet.com)  
[www.Detecto.com](http://www.Detecto.com)