

# EXAM ROOM PRODUCTS

From state-of-the-art sonar height rod scales with advanced electronics to traditional weighbeam physician scales



# DETECTO®

# DETECTO BRINGS CONVENIENCE AND ACCURACY TO YOUR EXAM ROOM WITH A VARIETY OF SCALES AND PRODUCTS YOU CAN RELY ON

Whether you're looking for the newest touchless scales that calculate BMI without pressing a single key, or a more traditional type weighbeam physician scale, DETECTO has a scale that's just right for your exam room. Our newest scales feature a touchless sonar height rod, and our traditional mechanical scale is still available and better than ever. Additional products include versatile mobile medical carts, waste receptacles, and steel glove box holders.

## APEX-SH

### Digital Physician Scale with Touchless Sonar Height Rod

DETECTO's APEX-SH series eye-level digital physician scales feature touchless sonar height rods and extra-wide, flat platforms measuring 17 in W x 17 in D x 2.75 in H. The scale automatically captures weight and height readings and displays both onscreen without pushing a single button. The scale features BMI calculation, 600 lb x 0.2 lb / 300 kg x 0.1 kg capacity, and large blue LCD weight readouts. Wi-Fi/Bluetooth models are also available.



## SOLO

### Eye-Level Digital Scale with Inline Mechanical Height Rod

The solo® eye-level clinical scale is the most economical, high-quality digital physician scale with height rod available. The scale features 550 lb x 0.2 lb / 250 kg x 0.1 kg capacity. The 14 inch x 15 inch platform is low-profile, and the inline mechanical height rod measures from 30 inches to 84.5 inches (76 cm to 214 cm). A USB port is also included.



## 339

### Traditional Weigh Beam Mechanical Scales with Height Rod, Handpost, and/or Wheels

DETECTO's world-renowned weighbeam physician scales feature a 10.5 in x 14.5 in (27cm x 37 cm) platform and dual-reading die-cast weighbeam which may be read from either side of the scale. Thirteen different models are available lb, lb/kg, and kg with capacities of 440 lb / 200 kg, 450 lb, or 200 kg, and included accessories such as handposts, height rods, and wheels.



## RC33669WHT

### Customizable RESCUE Steel Medical Carts

DETECTO's Rescue medical carts are sturdy and reliable... built to make sure you have the supplies you need where and when you need them. Our medical carts come in several colors and have your choice of five or six drawers that are ideal for storing a wide variety of supplies and can be secured. The ABS countertop is easy to clean and 5-inch wheels feature two parking locks and one steering lock. Optional accessories include IV poles, cardiac boards, catheter holders, drawer dividers, glove dispensers, sharps container holders and waste bins.



## GH3

### Universal Wall-Mount Steel Glove Box Holder

Conserve counter space with DETECTO's universal glove box holders. The open design allows you to see what size or type of glove is being dispensed. The holders accommodate most glove boxes and are universal enough to be used for tissue holders where counter or shelf space is limited. Models capable of holding one, two, or three boxes of gloves. Lifetime warranty.



## P-48R

### Step-On Cans, the Rugged Waste Receptacle

Named for their rugged toughness, Step-On Cans are perfect for laundry departments and receiving areas. Equipped with a heavy-duty foot pedal and linkage and a bag securing mechanism, they come in a variety of five different sizes and capacities (4, 6, 8, 12 and 25 gallon) and two different materials: stainless steel or a baked epoxy finish (available in red or white). The cans are flame retardant, have a self-closing lid for sanitation, and feature unitized construction.



*DETECTO reserves the right to improve, enhance, or modify features and specifications without prior notice.*

**SOLD BY:**

# DETECTO®

102 East Daugherty St., Webb City, MO 64870 USA  
Toll-Free: (800) 641-2008 E-mail: [detecto@cardet.com](mailto:detecto@cardet.com)  
[www.Detecto.com](http://www.Detecto.com)