

PH: (417) 673-4631 FAX: (417) 673-5001 WWW.DETECTO.COM



Ideal for home healthcare, uses for these lightweight, portable scales are almost limitless. With remote indicator and built-in carrying handle, the scale can go with you anywhere weighing is needed. The large, easy-to-clean platforms (2" H x 18" W x 14" D / 6 cm H x 46 cm W x 36 cm D) feature a black slip-proof mat and leveling eye. **LEGAL FOR TRADE.**



Certificate No. MV1 Display: 13-011 GP-400 Base: 94-098

MODEL	CAPACITY	SHIPPING WT.
GP-600-MV1	600 lb x 0.2 lb / 270 kg x 0.1 kg	31 lb / 14 kg
MV1PWR	AC Adapter	1 lb / 453 g



The scale's weight indicator includes an RS232 serial port, USB, or optional Wi-Fi Ethernet output for data transfer to a printer or a PC-based EMR (electronic medical record) software program. With the time/date feature, the scale can send the weight, time, date, BMI (Body Mass Index), and height directly to a printer or a PC.