


# GP-600-MV1

## Portable Low-Profile Healthcare Scale



Certificates:  
13-011 • 23-091



Ideal for home healthcare uses for these lightweight, portable scales are almost limitless. With remote indicator and built-in carrying handle, the scale can go with you anywhere weighing is needed. The large, easy-to-clean platforms (2" H x 18" W x 14" D) feature a black slip-proof mat and leveling eye. **LEGAL FOR TRADE.**

| MODEL             | CAPACITY                             | SHIPPING WEIGHT |
|-------------------|--------------------------------------|-----------------|
| <b>GP-600-MV1</b> | 600 lb x 0.2 lb /<br>270 kg x 0.1 kg | 31 lb / 14 kg   |
| <b>MV1PWR</b>     | AC Adapter                           | 1 lb / 453 g    |



The scale's weight indicator includes an RS232 serial port, USB, or optional Wi-Fi/Ethernet output for data transfer to a printer or a PC-based EMR (electronic medical record) software program. With the time/date feature, the scale can send the weight, time, date, BMI (Body Mass Index), and height directly to a printer or a PC.

# DETECTO®