


# HOME HEALTHCARE SCALES

Personal scales for bathrooms, mobile clinics, traveling nurses, telemedicine and general purpose weighing



**DETECTO**<sup>®</sup>

# HOME HEALTHCARE WEIGHING WITH A WIDE VARIETY OF PLATFORM SIZES AND FEATURES

DETECTO's home healthcare scales come with sturdy handles for easy transport, rubber mats that offer user stability, and a variety of displays to use that can be dependent on the patient's needs. The platforms are remarkably thin, allowing patients and home healthcare users to easily step onto the scale. The scales can weigh patients all the way up to 660 pounds, so most patients can be easily accommodated with most of the available home healthcare models. A variety of features and functions are available, including talking scales which offer audible confirmation of weights in multiple languages.

## SlimTALK

Talking Scale in English and Spanish with 400 lb / 180 kg Capacity



DETECTO's SlimTALK talking bathroom scales are ideal for the elderly, obese, blind, and telemedicine use. The scale is easy to use for the vision impaired and provides a clear, pleasant voice announcement of the weight in two languages: English or Spanish or OFF (switch selectable). The very low-profile platform is only 1.4 in/36 mm high for easy access. The bilingual SlimTALK scale features 400 lb x 0.1 lb / 180 kg x 0.1 kg capacity, 3 AAA batteries included, and large 1.5 in/38 mm high LCD digits for clear readouts.

## SlimTALK XL

Talking Scale in English and Spanish with 550 lb / 250 kg Capacity



The bilingual SlimTALK XL scale features 550 lb x 0.1 lb / 250 kg x 0.1 kg capacity, 3 AAA batteries included, large 1.5 in/38 mm high LCD digits for clear readouts, and 14.3 in W x 12.3 in D x 1.4 in H / 363 mm W x 312 mm D x 36 mm H overall dimensions. The step-on/off function allows for weighing without stooping to switch on the scales and conserves battery life. The SlimTALK XL has a durable plastic platform with ribbed black mat for weighing traction and user safety.

## SlimPRO

Premium Quality Low-Profile Scale That Speaks Two Different Languages



DETECTO's SlimPRO scale features 600 lb/272 kg capacity, Mother/Baby weighing mode, ribbed rubber mat for traction, backlit blue LCD, and an integral carrying handle. The SlimPRO's spacious platform size is 17 in W x 14.5 in D x 1.4 in H / 43 cm W x 37 cm D x 3 cm H. The scale is easy to use for the vision impaired and provides a clear, pleasant voice announcement of the weight in your choice of two different languages: English, Spanish, or OFF.

## D1130

Bathroom Dial Scales



Well-designed and quick-and-easy-to-read, DETECTO's ProHealth series scales are well-suited for locker rooms, hotels, fitness clubs, pharmacies, nutritional counseling centers, offices, and private households. DETECTO's D1130 bathroom scale features a large 7 in/18 cm dial, four colorful weight tracker markers, spacious platform that accommodates large shoe sizes (11.3 in x 11.8 in / 29 cm x 30 cm), non-skid rubber mat, 300 lb x 1 lb capacity, and adjusting scroll to set the zero point and balance.

## DR550C


Digital Portable Scale with Remote Readout



DETECTO's DR550C compact scale features 550 lb x 0.2 lb capacity, 12.5 in W x 11.9 in D x 2.2 in H / 318 mm W x 303 mm D x 58 mm H platform, simple four-key operation, slip-resistant texturized black rubber mat for traction and patient stability, and the coiled cable from scale base allows versatile positioning of the indicator, up to 90 inches away. The DR550C come with a UL listed 9V 300mA DC adapter (included) or may be powered by 4 AA batteries (not included) for mobile use. Optional black carrying case is available.

## DR660

Portable Bariatric Scale with Carrying Handle



DETECTO's DR660 portable bariatric scale features a lightweight spacious 21.5 in W x 17.2 in D x 1.97 in H (56 cm W x 44 cm D x 5 cm H cm) weighing platform, 660 lb x 0.5 lb / 300 kg x 0.2 kg capacity, 1 in/25 mm high LCD digits, RS232 serial output, integral carrying handle, and built-in rechargeable battery pack (7.2 VDC / 1200mA) which is able to provide eight hours of continuous use or 120 VAC 60Hz plug-in UL/CUL listed power supply (included). The display cable flexes from 34 in to 126 in (86 cm to 320 cm).

## GP-600-MV1

NTEP Legal for Trade for Athletic Wrestling or Dietary Clinic Weigh-in



Featuring battery power and a built-in carrying handle with a compact base design, these scales go anywhere. A large, easy-to-clean platform measuring 18 in D x 14 in W x 2 in H (46 cm D x 36 cm W x 6 cm H) features a black non-skid mat and bubble level for finding a level surface. The Lock/Release key can hold weights after stepping off the scale. The GP-600-MV1 is NTEP legal for trade and calculates Body Mass Index. Standard connectivity includes USB and RS232 serial, with additional models available that offer Wi-Fi/Bluetooth.

### Learn More:

<https://detecto.com/product/product-category/Home-Healthcare-Scales>

DETECTO reserves the right to improve, enhance, or modify features and specifications without prior notice.

### SOLD BY:



102 East Daugherty St., Webb City, MO 64870 USA  
Toll-Free: (800) 641-2008 E-mail: [detecto@cardet.com](mailto:detecto@cardet.com)  
[www.Detecto.com](http://www.Detecto.com)