

# ICU PRODUCTS

Medical carts, waste receptacles, and scales built to meet the demands of an ICU



**DETECTO**<sup>®</sup>

# THE STAKES ARE HIGH IN INTENSIVE CARE UNITS

To succeed you need equipment you can depend on for optimal patient satisfaction, like DETECTO's reliable and flexible medical carts, accurate scales, and tough waste receptacles.

Medical equipment from DETECTO is built tough to withstand the serious demands of an Intensive Care Unit. DETECTO medical carts are solid, secure, and customizable to perfectly fit your needs. Our hydraulic in-bed scales are battery or AC powered, with features that make it easy to provide accurate patient measurement for bed-bound patients. We also offer waste receptacles in many sizes and glove box holders built to last a lifetime.



## RC33669BLU-L RC33669RED-L

RESCUE Steel  
Medical Carts with  
Customizable  
Accessories

DETECTO's Rescue medical carts are sturdy and reliable... built to make sure you have the supplies you need where and when you need them. They come in several colors, with red selected for ER use and blue for Anesthesiology. Carts have your choice of five or six drawers that can be secured with a key or EMG breakaway lock. The ABS countertop is easy to clean, and 5-inch wheels feature two parking locks and one steering lock. Optional accessories include IV poles, cardiac boards, catheter holders, drawer dividers, glove dispensers, sharps container holders and waste bins.

## P-100R

Step-On Cans, the  
Rugged Waste Receptacle



Named for their rugged toughness, Step-On Cans are perfect for hospitals, doctors' offices, patients' rooms, and laboratory settings. Equipped with a heavy-duty foot pedal and linkage and a bag securing mechanism, they come in a variety of five different sizes and capacities (4, 6, 8, 12 and 25 gallon) and two different materials: stainless steel or a baked epoxy finish (available in red or white). The cans are flame retardant, have a self-closing lid for sanitation, and feature unitized construction.

## ID-7248SH-855RMP

Solace In-Floor Flush-Mounted  
Dialysis Scale



Solace series in-floor dialysis scales are flush-mounted platform scales that feature stainless steel decks and include a pit frame, trim ring, and interface cable to the included multi-color, touchscreen LCD indicator and printer. The flush surface eliminates tripping hazards and hallway accessibility issues. The Solace series is available with or without a handrail and all models include a recessed-mount indicator and ticket printer combination. Special dialysis weighing software is included with the scale. These in-floor scales with high 1,000 lb x 0.2 lb/450 kg x 0.1 kg capacity and a variety of sizes from 3 x 3 to 6 x 4 feet.

## IBFL500

Hydraulic Lift In-bed  
Scale for Patients  
up to 500 Pounds



DETECTO's in-bed scales feature flame-retardant, antibacterial 6 ft (1.8 m) long stretchers which are placed underneath the patient and hydraulically lifted inches from the bed for weighing, making them ideal for ICUs. They feature locking casters, an adjustable weighing arm, and steel construction. The ZERO key allows staff to easily subtract the weight of bedding for accurate patient measurement. The IBFL500 is AC or battery powered, EMR/EHR ready, and has a capacity of up to 500 lbs.

## IB800

Hydraulic Lift In-bed  
Scale for Patients  
up to 800 Pounds



DETECTO's IB800 scale is built with a base widening stance for added stability to handle a capacity of up to 800 pounds. It has a flame-retardant, antibacterial 6 ft (1.8 m) long stretcher which is placed underneath the patient and hydraulically lifted inches from the bed for weighing. The scale features locking casters, adjustable weighing arm, and steel construction. The ZERO key allows staff to easily subtract the weight of bedding for accurate patient measurement. AC or battery powered and EMR/EHR ready.

## GH3

Universal  
Wall-Mount  
Steel Glove Box  
Holders



Conserve counter space with DETECTO's universal glove box holders. The open design allows you to see what size or type of glove is being dispensed. The holders accommodate most glove boxes and are universal enough to be used for tissue holders where counter or shelf space is limited. Models capable of holding one, two, or three boxes of gloves. Lifetime warranty.

DETECTO reserves the right to improve, enhance,  
or modify features and specifications without prior notice.

**SOLD BY:**

# DETECTO®

102 East Daugherty St., Webb City, MO 64870 USA  
Toll-Free: (800) 641-2008 E-mail: [detecto@cardet.com](mailto:detecto@cardet.com)  
[www.Detecto.com](http://www.Detecto.com)