



model **MLM** **INFANT STADIOMETER**

MECHANICAL LENGTH MEASURING



Detecto
Model MLM
Mechanical Infant
Stadiometer



- High-contrast, easy-to-read measurements
- Graduations in 0.125 in / 1 mm increments
- Measurement ranging from 2.5 to 39 in and 65 to 1,000 mm
- May be used with babies and infants up to 2 years old
- Durable construction which requires no electricity to operate


Pediatricians' Offices • Clinics • Hospitals

DETECTO's mechanical infant stadiometer is a durable length measuring device which ensures years of trouble-free operation. The model MLM's portable slider features soft ergonomic edges to protect young patients from unnecessary injury and maximize patient comfort while measuring.

model **MLM** INFANT STADIOMETER



Lightweight and Easy to Use



Model	MLM
Measuring Range	2.5 - 39 inches / 65 - 1,000 mm
Graduation	0.125 inch / 1 mm
Accuracy	± 0.125 inch / ± 1 mm
UPC Code	809161188301
Shipping Weight	2 lb / .9 kg

The instrument's smart design

allows the measurement to be retained after usage freeing the medical professional to attend to the patient first before recording the measurement. The light cream color is neutral to compliment all medical settings and the high-contrast measurements are easy to read.



Easy pushbutton extension of the measurement caliper



The calipers of the portable infant stadiometer conveniently fold in for easy transport or storage.

Detecto Scale reserves the right to improve, enhance, or modify features and specifications without prior notice.

 **DETECTO**[®]
A Division of Cardinal Scale Manufacturing Co.
203 E. Daugherty, Webb City, MO 64870 USA
Ph: 417-673-4631 or 1-800-641-2008 Fax: 417-673-5001
www.Detecto.com

Sold By: