

# NICU PRODUCTS

Scales, medical carts, and equipment to help you  
care for your most fragile patients



# NICU'S DELICATE PATIENTS REQUIRE THE UTMOST SUPPORT AND PRECISE EQUIPMENT

**DETECTO's versatile and accurate scales get you the information you need quickly and easily for recording patient measurements.**

For more than 120 years, DETECTO has been making medical scales to serve the needs of doctors, nurses and healthcare professionals around the world. Our diverse line of hospital scales will meet the demands of the NICU department, where speed and accuracy are necessary to care for the most fragile patients. DETECTO's mobile baby scale carts and waste receptacles are built to perfectly serve your needs and perform reliably for a lifetime.

## DLM

**Digital Length Measuring Device**



DETECTO's DLM Digital Length Measuring Device may be used stand-alone in the NICU or with the model MB130 pediatric scale. The lightweight design and battery power make it ideal for easy transport from one patient room to another. Its bold LCD is easy to read, and it features convenient unit switching, a hold feature, and overall dimensions of 11.4 W x 2.7 H x 24.4 L inch (29 W x 7 H x 62 L cm). Powered by two 1.5V AAA batteries (included). Length range: 13.75 to 31.5 inch (35 to 80 cm).

## 6745

**Digital Baby Scale, Weighs up to 30 lb**



DETECTO's 6745 baby scale was designed for use in the NICU, nursery, delivery room and pediatric practices where high accuracy and dependability is paramount for weighing the smallest of patients. Step-by-step instructions are displayed right on the keypad for ease of use. The baby scale features metric lockout, integral measuring tape (up to 22 inch/56 cm), 30 lb x 0.1 oz / 15 kg x 0.005 kg capacity, 1 inch (25 mm) high LCD, and RS232 serial output for connectivity with a printer or PC. The spacious weighing tray measures 26 in W x 12 in D / 66 cm W x 31 cm D.

## MB130

**Digital Pediatric Scale with Extra-Wide Weighing Tray**



This flexible scale lets you weigh patients from neonatal up to toddler with its 44lb/20kg capacity. The extra-large weighing tray keeps babies safe and comfortable, and the 5 gram accuracy is great for weighing your smallest patients. A measuring tape is built in, with a digital length measuring rod available. The Breast Milk Intake function compares the weight of a baby before and after feeding, and then alternates the display between the baby's current weight and the change from the first weight, indicating the amount of milk intake of the baby.

## WPS12UT

**Digital Diaper Scale is Submersible Underwater**



The WPS12UT digital scale with utility bowl is IP67-rated for water protection and can be completely submersed underwater for easy cleaning. This stainless steel scale is ideal for weighing soiled or wet diapers, chux pads, lap sponges, and other materials. The 4.6 in / 11.7 cm deep bowl aids sample retention and allows easy clean-up. The scale's digital speed and accuracy means you'll be able to view weights instantly with up to a fine 0.1 ounce or 1 gram division size.

## P-32R

**Step-On Cans, the Rugged Waste Receptacle**



Named for their rugged toughness, Step-On Cans are perfect for nearly anywhere, with the white 32-quart unit just right for Neonatal departments. Equipped with a heavy-duty foot pedal and linkage and a bag securing mechanism, they come in a variety of five different sizes and capacities (4, 6, 8, 12 and 25 gallon) and two different materials: stainless steel, or a baked epoxy finish (available in red or white). The cans are flame retardant, have a self-closing lid for sanitation, and feature unitized construction.

## AP-4KD

**Compact Wet Diaper and Lap Sponge Scale**



This compact scale is perfect for measuring fluid output, and the removable stainless steel tray is easy to clean. The instantaneous display of .5 inch high red digits makes weighing quick and easy. The handy push-button tare makes things simple. The 4,000 g capacity weighs accurately down to 1 gram.

## RC33669WHT

**White RESCUE Baby Scale Carts**



DETECTO's Rescue series baby scale perfectly accommodates many of our baby scales and length measuring rods for a mobile NICU measurement station that can go anywhere. The steel rolling cart with a white powder paint features a push handle and two parking locks plus one steering lock for optimal patient safety. The top shelf is made of ABS molded plastic and will securely grip the baby scale's four rubber feet without sliding for equipment stability and patient safety.

*DETECTO reserves the right to improve, enhance, or modify features and specifications without prior notice.*



102 E. Daugherty, Webb City, MO 64870 USA  
Ph: 417-673-4631 or 1-800-641-2008 Fax: 417-673-2153  
[www.Detecto.com](http://www.Detecto.com)

**SOLD BY:**