

RESCUE
medical carts



PEDIATRIC CART

9 Multi-Color Drawers



*DETECTO model
RC33333336PED*

DETECTO's Rescue series Pediatric cart is designed for use with the proven clinical methodology for treating life-threatening emergencies in pediatric patients. When seconds count, the color-coded drawers enable doctors and clinicians to quickly identify critical supplies while responding to pediatric emergencies.

Contents are secured with individual plastic seals for each specific drawer or a keyed lock which provides immediate access to all drawers and emergency supplies.



DETECTO



DETECTO Rescue Series Pediatric Medical Cart with Push Handle and User-Friendly Design

DETECTO's Rescue series pediatric carts come in 9 multi-colored drawers using the Broselow® color-coding method and have 5-inch diameter wheels which feature two parking locks and one steering lock. Rescue carts are manufactured from quality steel cart body construction that is powder-coat painted for durability with ABS countertops that are easy to clean. ABS bumpers protect all four corners of the medical carts. The push handle offers a user-friendly design.



- White frame with 9 multicolor drawers
- Drawer Sizes: 3 in cream / 3 in pink/red / 3 in lavender / 3 in yellow / 3 in white / 3 in light blue / 3 in bright orange / 3 in fresh green / 6 in cream
- Each drawer has an individual breakaway lock
- Keyed lock for all drawers
- Numerous optional accessories available
- Push handle with user-friendly design
- Auto-return rail on drawers
- ABS countertop surface is easy to clean
- Quality steel cart body construction
- 2 parking locks and 1 steering lock on casters
- 5-in diameter wheels
- ABS bumpers on each of the 4 corners



Individual drawer locks and universal keyed lock for all drawers. Optional CAPS red plastic seals are available that come in qty. 50, individually numbered.

SPECIFICATIONS	
Model	RC33333336PED
Overall Dimensions (without push handle)	26.9 in W x 21.5 in D x 45 in H / 68.3 cm W x 54.6 cm D x 114.3 cm H
Overall Dimensions (with push handle)	29.8 in W x 21.5 in D x 47 in H / 75.7 cm W x 54.6 cm D x 119.4 cm H
Inner Drawer Dimensions (3-inch-high drawers)	22.3 in W x 15.5 in D x 2.5 in H / 56.6 cm W x 39.4 cm D x 6.4 cm H
Inner Drawer Dimensions (6-inch-high drawers)	22.3 in W x 15.5 in D x 5.5 in H / 56.6 cm W x 39.4 cm D x 14.0 cm H
Shipping Weight	161 lb / 73 kg
Shipping Dimensions	28 in x 22 in x 46 in / 71 cm x 56 cm x 117 cm
Options	CAPS red plastic seals (qty. 50, individually wrapped) and a wide range of optional accessories including bin organizers, cardiac board, defibrillator shelf, IV pole, drawer dividers, sharps container holder, catheter holder, oxygen tank holder, and waste bin.



102 East Daugherty St., Webb City, MO 64870 USA
Toll-Free: (800) 641-2008 E-mail: detecto@cardet.com
www.Detecto.com

SOLD BY: