




# ProHealth Dial Scales

D1130 / D1130K



**PERSONAL BATHROOM SCALES**

**PROHEALTH**

*Well-designed and quick-and-easy-to-read, Detecto's ProHealth Series Scales are well-suited for locker rooms, hotels, fitness clubs, pharmacies, nutritional counseling centers, offices, and private households. They feature fast weighing with the weight display immediately registering on the large 7-inch dial after stepping on the scale's slip-resistant platform.*

# PROHEALTH SCALES



- Both Pounds and Kilograms Versions are Available.
- Easy-to-Read, Large 7-Inch Dial Registers Weight Quickly.
- Features Adjusting Scroll to Set the Zeropoint and Balance.
- Spacious Platform Accommodates Large Shoe Sizes.
- Low-Profile, Compact Design Requires Little Floor Space.

## DETECTO MODEL D1130 / D1130K PROHEALTH SCALES

A slim, low-profile design, economical pricing, and non-skid vinyl mat are just a few features of the Detecto D1130 and D1130K scales - proven performers for budget-minded customers. They feature a white, baked-enamel finish and personal, colorful indicators to mark one

person's or a whole family's weights. The large 7-inch-diameter dial allows for easy readability even for tall people. Designed for personal use, fitness centers, and health care facilities. Both pounds and kilograms versions are available.



model D1130

### D1130 SPECIFICATIONS

<b>CAPACITY</b>	300 lb x 1 lb
<b>DIMENSIONS (OVERALL):</b>	12" W x 16.75" D x 2.5" H 305 mm W x 425 mm D x 64 mm H
<b>DIMENSIONS (PLATFORM):</b>	11.25" W x 11.75" D 286 mm W x 298 mm D
<b>DIAL DIAMETER:</b>	7" / 178 mm
<b>SHIPPING WEIGHT:</b>	9 lb / 4 kg



model D1130K

### D1130K SPECIFICATIONS

<b>CAPACITY</b>	130 kg x 500 g
<b>DIMENSIONS (OVERALL):</b>	12" W x 16.75" D x 2.5" H 305 mm W x 425 mm D x 64 mm H
<b>DIMENSIONS (PLATFORM):</b>	11.25" W x 11.75" D 286 mm W x 298 mm D
<b>DIAL DIAMETER:</b>	7" / 178 mm
<b>SHIPPING WEIGHT:</b>	9 lb / 4 kg

*Detecto Scale reserves the right to improve, enhance, or modify features and specifications without prior notice.*

 **DETECTO**<sup>®</sup>  
A Division of Cardinal Scale Manufacturing Co.  
102 E. Daugherty, Webb City, MO 64870 USA  
Pb: 417-673-4631 or 1-800-641-2008 Fax: 417-673-2153  
[www.Detecto.com](http://www.Detecto.com)

**Sold By:**