



360-DEGREE ROTATION FOR TIGHT CORNERS

RESCUE medical carts



- Auto-return rail on drawers
- Keyed or EMG breakaway central lock secures all drawers at the same time
- Push handle with userfriendly design
- Colors available for ER (red), Anesthesiology (blue), Isolation (yellow), and General Purpose (white) departments
- ABS countertop surface is easy to clean
- Quality aluminum cart body construction
- ABS bumpers on each of the 4 corners
- · 5-in diameter wheels
- Numerous optional accessories available such as IV poles, cardiac boards, catheter holders, drawer dividers, glove dispensers, sharps container holders, and waste bins
- Drawer dimensions: 22.4 in W x 15.6 in D / 56.9 cm W x 39.6 cm D
- Cart dimensions: 27.8 in W x 21.3 in D x 41.8 in H (without guide handle) / 70.6 cm W x 54.1 cm D x 106.2 cm H



Hygienic Wipe-down Surfaces

The Rescue series carts are built with the busy clinical user in mind and feature an ABS countertop with smooth molded contours that are easy to clean for hygienic wipe-down disinfection. A central lock secures all the drawers at the same time for convenience.



Versatile Wheel Locks

The Rescue carts feature 2 parking locks and 1 directional steering lock on the wheels. The 5-inch-diameter wheels are built to last under heavy-duty clinical daily use. ABS bumpers protect all four corners of the medical carts.



Multiple Accessories Available

The Rescue series features a number of accessories available (like these drawer dividers) or you may purchase a fully-loaded cart pre-stocked with typical accessories needed for ER or Anesthesiology departments.

EMERGENCY ROOM CARTS ANESTHESIOLOGY CARTS 5-Drawer Cart 5-Drawer Cart Model: CAIP Model: RC33669BLU **Description:** Chrome IV Pole Model: RC33669RED Description: White **Description:** White Dimensions: Frame with (5) Red Frame with (5) Blue 3"/3"/6"/9" Drawers 3"/3"/6"/9" Drawers with EMG Lock with Keyed Lock Shipping Wt: 132 lb Shipping Wt: 132 lb Model: CADS 5-Drawer Cart 5-Drawer Cart Model: RC33669RED-L Model: RC33669BLU-L Dimensions: Description: White **Description:** White Frame Frame Loaded Cart with Loaded Cart with (5) Blue 3"/3"/6"/9" Drawers with (5) Red 3"/3"/6"/6"/9" Keyed Lock Drawers with EMG Lock Included Accessories: **Included Accessories:** IV pole, defibrillator shelf, sharps container holder, cardiac board, oxygen 12-bin organizer with Shipping Wt: 6 lb tank holder, sharps accessory bridge, Model: CACB container holder, and catheter holder, waste bin, and two two accessory rails accessory rails Dimensions: Shipping Wt: Cart: 132 lb Shipping Wt: Cart: 132 lb Accessory Box: 22 lb Accessory Box: 18 lb 6-Drawer Cart 6-Drawer Cart Shipping Wt: 7 lb Model: RC333369RED Model: RC333369BLU Description: White **Description:** White Frame with (6) Red Frame with (6) Blue 3"/3"/3"/6"/9" Drawers 3"/3"/3"/6"/9" Drawers Dimensions: with EMG Lock with Keyed Lock Shipping Wt: 132 lb Shipping Wt: 132 lb Shipping Wt: 4 lb 6-Drawer Cart 6-Drawer Cart Model: RC333369RED-L Model: RC333369BLU-L **Description:** White Description: White Model: CARCDSC **Description:** Sharps Container Frame Loaded Cart with Frame with (6) Blue (6) Red 3"/3"/3"/3"/6"/9" 3"/3"/3"/6"/9" Drawers Drawers with EMG Lock with Keyed Lock Dimensions: Included Accessories: **Included Accessories:** IV pole, defibrillator shelf, sharps container holder, cardiac board, oxygen 12-bin organizer with Shipping Wt: 2 lb tank holder, sharps accessory bridge, container holder, and catheter holder, waste bin, and two two accessory rails Model: CAB12 accessory rails Shipping Wt: Cart: 132 lb Shipping Wt: Cart: 132 lb Accessory Box: 22 lb Accessory Box: 18 lb Dimensions: GENERAL PURPOSE / UTILITY CARTS **ISOLATION MEDICAL CARTS** 5-Drawer Cart 5-Drawer Cart Model: RC33669WHT Model: RC33669YEL **Description:** White **Description:** White Frame with (5) White Frame with (5) Yellow Model: CARCCH 3"/3"/6"/6"/9" Drawers 3"/3"/6"/6"/9" Drawers with Keyed Lock with Keyed Lock Accessory Rail Shipping Wt: 132 lb Shipping Wt: 132 lb Dimensions: Shipping Wt: 2 lb 6-Drawer Cart 6-Drawer Cart Model: RC333369WHT Model: RC333369YEL Model: CARCWB



Description: White Frame with (6) White 3"/3"/3"/6"/9" Drawers with Keyed Lock

Shipping Wt: 132 lb



Description: White Frame with (6) Yellow 3"/3"/3"/6"/9" Drawers with Keyed Lock


Shipping Wt: 132 lb




12.1 in W x 1.5 in D x 27.6 in H (46.3 in H fully extended) / 30.7 cm W x 3.8 cm D x 70.1 cm H (117.6 cm H fully extended)

Shipping Wt: 4 lb

Description: Aluminum Defibrillator Shelf





17.3 in W x 0.3 in D x 20.1 in H / 43.9 cm W x 0.8 cm D x 51.1 cm H

Model: CARCOH



14.7 cm W x 12.4 cm D x 40.4 cm H



Holder with Accessory Rail

11.0 in W x 5.5 in D x 6.4 in H / 27.9 cm W x 14.0 cm D x 16.3 cm H

ANESTHESIOLOGY LOADED CART ACCESSORIES

Description: 12-Bin Organizer with Accessory Bridge

23.6 in W x 3.3 in D x 9.8 in H / 60 cm W x 8.4 cm D x 24.9 cm H

Shipping Wt: 6 lb PLUS 2 lb for accessory bridge

Description: Catheter Holder with

4.0 in W x 4.6 in D x 18.9 in H / 10.2 cm W x 11.7 cm D x 48 cm H



Description: Waste Bin with Accessory Rail

Dimensions:

10.4 in W x 7.0 in D x 17.3 in H / 26.4 cm W x 17.8 cm D x 43.9 cm H Shipping Wt: 4 lb

Model: CAB05 Organizer with Accessory Bridge Dimensions:

Description: 5-Bin

23.5 in W x 4.8 in D x 5.9 in H / 59.7 cm W x 12.2 cm D x 15 cm H

Shipping Wt: 5 lb PLUS 2 lb for accessory bridge



Model: CARCGD

ADDITIONAL RESCUE CART ACCESSORIES

Description: Glove Dispenser with Accessory Rail

Dimensions: 10.5 in W x 4.5 in D x 5.3 in H / 26.7 cm W x 11.4 cm D x 13.5 cm H Shipping Wt: 2 lb

Model: CARCMB

Description: Multi-Storage Bin with Accessory Rail

Dimensions:

13.1 in W x 5.6 in D x 7.1 in H / 33.3 cm W x 14.2 cm D x 18 cm H

Shipping Wt: 3 lb



Model: CAPS Description: 50 Red Plastic Seals, Individually Numbered

Shipping Wt: 1 lb



Model: CARCDS3

Description: Cart Divider Set for 3" Drawer, Includes 3 short and 3 long dividers for 16 compartment capability

Shipping Wt: 2 lb



Model: CARCDS6 **Description:** Cart Divider

Set for 6" Drawer, Includes 2 short and 2 long dividers for 9 compartment capability

Shipping Wt: 2 lb

CUSTOM CARTS AND PEDIATRIC CARTS:

If you don't see a standard medical cart model in this brochure that suits your application, DETECTO can custom build it for your clinical facility. One of our representatives will be happy to discuss your particular needs, and we can provide medical carts built exactly to your specifications.





Detecto Scale reserves the right to improve, enhance, or modify features and specifications without prior notice.

DETECTO www.Detecto.com

102 E. Daugherty, Webb City, MO 64870 USA Ph: 417-673-4631 or 1-800-641-2008 Fax: 417-673-2153 www.Detecto.com

SOLD BY: