



# 360-DEGREE ROTATION FOR TIGHT CORNERS

DETECTO

# RESCUE medical carts

- Auto-return rail on drawers
- Keyed or EMG breakaway central lock secures all drawers at the same time
- Push handle with userfriendly design
- Colors available for ER (red), Anesthesiology (blue), Isolation (yellow), and General Purpose (white) departments
- ABS countertop surface is easy to clean
- Quality steel cart body construction
- ABS bumpers on each of the 4 corners
- 5-in diameter wheels
- Numerous optional accessories available such as IV poles, cardiac boards, catheter holders, drawer dividers, glove dispensers, sharps container holders, and waste bins
- Drawer dimensions: 22.4
  in W x 15.6 in D / 56.9 cm
  W x 39.6 cm D
- Cart dimensions: 27.8 in W x 21.3 in D x 41.8 in H (without guide handle) / 70.6 cm W x 54.1 cm D x 106.2 cm H



Model RC33669BLU-L

## **Hygienic Wipe-down Surfaces**

The Rescue series carts are built with the busy clinical user in mind and feature an ABS countertop with smooth molded contours that are easy to clean for hygienic wipe-down disinfection. A central lock secures all the drawers at the same time for convenience.



•

## Versatile Wheel Locks

The Rescue carts feature 2 parking locks and 1 directional steering lock on the wheels. The 5-inch-diameter wheels are built to last under heavy-duty clinical daily use. ABS bumpers protect all four corners of the medical carts.



360

## **Multiple Accessories Available**

The Rescue series features a number of accessories available (like these drawer dividers) or you may purchase a fully-loaded cart pre-stocked with typical accessories needed for ER or Anesthesiology departments.

#### **EMERGENCY ROOM CARTS**

5-Drawer Cart

5-Drawer Cart

Model: RC33669RED-L

Frame Loaded Cart with

Drawers with EMG Lock

IV pole, defibrillator shelf,

**Included Accessories:** 

cardiac board, oxygen

tank holder, sharps

with EMG Lock

6-Drawer Cart

Shipping Wt: 132 lb

Description: White

(5) Red 3"/3"/6"/6"/9"



#### Model: RC33669RED Description: White Frame with (5) Red 3"/3"/6"/6"/9" Drawers with EMG Lock Shipping Wt: 132 lb



## ANESTHESIOLOGY CARTS

5-Drawer Cart Model: RC33669BLU Description: White Frame with (5) Blue 3"/3"/6"/6"/9" Drawers with Keyed Lock Shipping Wt: 132 lb

## 5-Drawer Cart

Model: RC33669BLU-L Description: White Frame Loaded Cart with (5) Blue 3"/3"/6"/6"/9" Drawers with Keyed Lock Included Accessories:

sharps container holder, 12-bin organizer with accessory bridge, catheter holder, waste bin, and two accessory rails Shipping Wt: Cart: 132 lb Accessory Box: 18 lb

6-Drawer Cart Model: RC333369BLU Description: White Frame with (6) Blue 3"/3"/3"/3"/6"/9" Drawers with Keyed Lock Shipping Wt: 132 lb

6-Drawer Cart Model: RC333369BLU-L Description: White Frame with (6) Blue 3"/3"/3"/3"/6"/9" Drawers with Keyed Lock **Included Accessories:** sharps container holder, 12-bin organizer with accessory bridge, catheter holder, waste bin, and two accessory rails Shipping Wt: Cart: 132 lb Accessory Box: 18 lb

### **ISOLATION MEDICAL CARTS**

5-Drawer Cart Model: RC33669YEL Description: White Frame with (5) Yellow 3"/3"/6"/6"/9" Drawers with Keyed Lock Shipping Wt: 132 lb

Model: RC333369YEL Description: White Frame with (6) Yellow 3"/3"/3"/3"/6"/9" Drawers with Keyed Lock Shipping Wt: 132 lb

#### **ER LOADED CART ACCESSORIES**

Model: CAIP Description: Chrome IV Pole Dimensions: 12.1 in W x 1.5 in D x 27.6 in H (46.3 in H fully extended) / 30.7 cm W x 3.8 cm D x 70.1 cm H (117.6 cm H fully extended)

Shipping Wt: 4 lb

Model: CADS Description: Aluminum Defibrillator Shelf

Dimensions: 13.1 in W (16.4 in W fully extended) x 13.3 in D x 16.1 in H (26 in H fully extended) / 33.3 cm W (41.7 cm W fully extended) x 33.8 cm D x 40.9 cm H (66 cm H fully extended) Shipping Wt: 6 lb

Model: CACB **Description:** Clear Acrylic Cardiac Board with Handle

Dimensions: 17.3 in W x 0.3 in D x 20.1 in H / 43.9 cm W x 0.8 cm D x 51.1 cm H Shipping Wt: 7 lb

Model: CARCOH Description: Oxygen Tank Holder with Accessory Rail Dimensions:

5.8 in W x 4.9 in D x 15.9 in H / 14.7 cm W x 12.4 cm D x 40.4 cm H Shipping Wt: 4 lb


#### ER AND ANESTHESIOLOGY LOADED CART ACCESSORY

Model: CARCDSC Description: Sharps Container Holder with Accessory Rail Dimensions: 11.0 in W x 5.5 in D x 6.4 in H / 27.9 cm W x 14.0 cm D x 16.3 cm H Shipping Wt: 2 lb

ANESTHESIOLOGY LOADED CART ACCESSORIES

Model: CAB12 Description: 12-Bin Organizer with Accessory Bridge

> Dimensions: 23.6 in W x 3.3 in D x 9.8 in H / 60 cm W x 8.4 cm D x 24.9 cm H Shipping Wt: 6 lb PLUS 2 lb for accessory bridge



Model: CARCWB Description: Waste Bin with Accessory Rail Dimensions: 10.4 in W x 7.0 in D x 17.3 in H / 26.4 cm W x 17.8 cm D x 43.9 cm H Shipping Wt: 4 lb



container holder, and two accessory rails Shipping Wt: Cart: 132 lb Accessory Box: 22 lb 6-Drawer Cart Model: RC333369RED Description: White Frame with (6) Red 3"/3"/3"/3"/6"/9" Drawers



Description: White Frame Loaded Cart with (6) Red 3"/3"/3"/3"/6"/9" Drawers with EMG Lock Included Accessories: IV pole, defibrillator shelf, cardiac board, oxygen tank holder, sharps container holder, and two accessory rails

Shipping Wt: Cart: 132 lb Accessory Box: 22 lb

#### **GENERAL PURPOSE / UTILITY CARTS**



5-Drawer Cart Model: RC33669WHT Description: White Frame with (5) White 3"/3"/6"/6"/9" Drawers with Keyed Lock Shipping Wt: 132 lb

6-Drawer Cart Model: RC333369WHT Description: White Frame with (6) White 3"/3"/3"/3"/6"/9" Drawers with Keyed Lock Shipping Wt: 132 lb



6-Drawer Cart



ADDITIONAL RESCUE CART ACCESSORIES					
	Model: CAB05 Description: 5-Bin Organizer with Accessory Bridge Dimensions: 23.5 in W x 4.8 in D x 5.9 in H / 59.7 cm W x 12.2 cm D x 15 cm H Shipping Wt: 5 lb PLUS 2 lb for accessory bridge		Model: CARCGD Description: Glove Dispenser with Accessory Rail Dimensions: 10.5 in W x 4.5 in D x 5.3 in H / 26.7 cm W x 11.4 cm D x 13.5 cm H Shipping Wt: 2 lb		Model: CARCMB Description: Multi-Storage Bin with Accessory Rail Dimensions: 13.1 in W x 5.6 in D x 7.1 in H / 33.3 cm W x 14.2 cm D x 18 cm H Shipping Wt: 3 lb
A \$ 05.572	Model: CAPS Description: 50 Red Plastic Seals, Individually Numbered Shipping Wt: 1 lb		Model: CARCDS3 Description: Cart Divider Set for 3" Drawer, Includes 3 short and 3 long dividers for 16 compartment capability Shipping Wt: 2 lb		Model: CARCDS6 Description: Cart Divider Set for 6" Drawer, Includes 2 short and 2 long dividers for 9 compartment capability Shipping Wt: 2 lb

DETECTO

8


#### CUSTOM CARTS AND PEDIATRIC CARTS:

If you don't see a standard medical cart model in this brochure that suits your application, DETECTO can custom build it for your clinical facility. One of our representatives will be happy to discuss your particular needs, and we can provide medical carts built exactly to your specifications.

> Detecto Scale reserves the right to improve, enhance, or modify features and specifications without prior notice.

.

SOLD BY:



**DETECTO** 

www.Detecto.com

© Copyright 2024 Cardinal Scale Mfg. Co. • Printed in USA • CAR/00/0124/C270A