

SOLO-RI

LOW-PROFILE SCALE WITH REMOTE INDICATOR



solo
DIGITAL CLINICAL SCALES®



- Low-profile scale base with 7.25 ft / 2.2 m cable to remote display
- Display may be wall-mounted or set on countertop (bracket included)
- Body Mass Index (BMI) calculation to determine health status
- 14.2 in W x 15.8 in D x 2.4 in H / 36.0 cm W x 40.2 cm D x 6.1 cm H platform size
- 550 lb x 0.2 lb / 250 kg x 0.1 kg capacity

MEDICAL-GRADE ACCURACY:

WEIGHT, HEIGHT, AND BODY MASS INDEX

SOLO-RI Portable Clinical Scales with Space-Saving Compact Design for Flexible Use

- 550 lb x 0.2 lb / 250 kg x 0.1 kg capacity
- Weight, Height, and Body Mass Index
- Low-profile platform (only 2.4 in / 6.1 cm high)
- 7.25 ft (2.2 m) cable from scale base to display to allow placing display wherever it's easiest to read
- Mother/Baby mode factors the weight of infants and toddlers held by an adult
- Wall or desk mount bracket for the display (included)
- Manual input height entry from 2 ft 0.0 in - 7 ft 10.5 in by 0.5 in (60.0 cm - 240.0 cm by 1cm)
- Audible beep confirmation (short beep when a key is pressed, and a long beep for when weight is stable and locked)
- Selectable Auto Shut-Off feature to conserve batteries (up to 9 minutes)
- Auto Weight Lock feature on stable weight
- 6 AA batteries included. SOLO-RI-AC model comes with AC adapter
- Included carpet feet improve weighing accuracy on carpeting

The SOLO-RI provides BMI (Body Mass Index) calculation by conveniently selecting your height while weighing with the easy-to-use scroll arrow keys.



DETECTO's SOLO-RI is portable and light weight for easily changing weighing locations based on patient needs.

SOLO
DIGITAL CLINICAL SCALES



The 7.25 ft/2.2 m cable from the scale base to the indicator allows you to set the display wherever it's most convenient.



The SOLO-RI's remote display may be positioned on countertops, desks, or mounted on the wall for flexible weighing use.

Countertop bracket with rubber feet included for level surfaces



DETECTO
www.Detecto.com

*The Name Brand Preferred
by Clinical Professionals
Around the World*



Stylish Indicator

The SOLO-RI indicator features 8 simple buttons that are easily identifiable in any language for use around the world.



OneWeigh Units Locking

OneWeigh units locking into Pounds or Kilograms avoids dosage errors based on incorrect units.



Wall-Mounting Bracket

A multifunctional bracket is included to mount the display on a wall, or the bracket (and included feet) can be used to place it on countertop.



Weight Hold Key

The SOLO-RI's Hold/Release key lets you hold a patient weight onscreen until you can record it in their recordkeeping.



EMR/EHR-Ready

The scale's USB port allows you to output measurement data to EMR/EHR software by pressing the Print key.



Wide, Flat Platform

The oversized 14 x 15 inch platform with ribbed rubber mat is only 2.4 inches high making it easy for patients to step on/off.



Everything Included

The SOLO-RI comes ready for weighing right out of the box with included batteries, manual, wall-mount and desktop bracket, screws, and carpet feet.



Battery Power

No power outlet available, no problem! The SOLO-RI comes with 6 AA batteries for weighing anywhere (model SOLO-RI-AC comes with AC adapter).



Model
SOLO

DETECTO's
SOLO-RI is also
available with an
eye-level column
as model SOLO



Model
SOLO-RI

Optional
carrying case
with internal
Velcro straps
available (model
PRODOC-CASE)



Model	SOLO-RI	SOLO-RI-AC
Capacity	550 lb x 0.2 lb / 250 kg x 0.1 kg	
Platform Size	14.2 in W x 15.8 in D x 2.4 in H / 36.0 cm W x 40.2 cm D x 6.1 cm H	
Remote Display Size (with mounting bracket)	7.8 in W x 5.1 in D x 3.0 in H / 19.8 cm W x 12.9 cm D x 7.7 cm H	
Weight Units	Pounds or Kilograms (selectable upon initialization and then locked in)	
Power	(6) AA size Alkaline batteries (included)	(6) AA size Alkaline batteries (included) and medical device 9V AC/DC wall plug-in adapter (DETECTO part number PD-ACKIT, includes UK plug, EU plug, and USA plug)
Display	Five-digit, seven-segment, 1.0-inch (25-mm) high LCD digits	
Operating Temperature	50 to 104 °F / +10 to +40 °C	
Function Keys	ON/OFF, ZERO, BMI HEIGHT, ENTER, LOCK/RELEASE, PRINT, UP ARROW/Mother and Baby, DOWN ARROW	
Annunciators	ZERO, LOCK, STA, HOLD, LB, KG, CM, FT/IN, HEIGHT, WEIGHT, BMI, AM, PM	
Interface	USB (ASCII Protocol)	
Height Units	Inches or centimeters	
Manual Input Height	2 ft 0.0 in – 7 ft 10.5 in by 0.5 in / 60.0 cm - 240.0 cm by 1 cm	
Minimum Weighing	11 lb / 5 kg	
Selectable Auto Shut-Off	Up to 9 minutes	
Construction	base – black ABS; mat – black PVC	
Product Weight	11 lb / 5 kg	12 lb / 5.5 kg
Shipping Weight	14 lb / 6.4 kg	15 lb / 6.9 kg
Package Dimensions	24 in W x 17 in D x 7 in H / 61 cm W x 43 cm D x 18 cm H	
UPC Code	809161330304	809161330700

Detecto Scale reserves the right to improve, enhance,
or modify features and specifications without prior notice.

SOLD BY:

 **DETECTO**
www.Detecto.com

102 E. Daugherty, Webb City, MO 64870 USA
Ph: 417-673-4631 or 1-800-641-2008 Fax: 417-673-2153
www.Detecto.com