



Slim PRO

# LOW-PROFILE HEALTHCARE SCALES



DETECTO's SlimPRO is a high-capacity, low-profile digital scale for clinical or personal weighing use. The spacious, extra-wide weighing platform easily accommodates patients up to 440 lb/200 kg.

## THIN, LOW-PROFILE HEALTHCARE WEIGHING

As the name implies, DETECTO's SlimPRO is an incredibly thin 1.4-in/3-cm high scale that contains loads of weighing performance and durability for clinical-grade usage weighing heavy patients.



## RS232 Serial Connectivity

**The SlimPRO goes a step above the competition** by providing a RS232 serial output port for connectivity with a printer or PC. This allows the SlimPRO to output weight data for efficient patient recordkeeping. DETECTO offers an optional cable for the serial output (shown left). The SlimPRO may be connected to a PC for ASCII format weight data to output automatically on a stable weight reading, or this data may be sent directly to a printer for a hard-copy receipt of the weight reading.



### Serial Interface Cable

DETECTO Part No. 6600-1940

- When connecting to a DB9 printer connection, a gender changer is required.
- When connecting to a PC, if the PC does not have a DB9 connection, then a USB/DB9 converter is required.



### P50 PRINTER

DETECTO's optional model P50 direct thermal printer may be connected to the SlimPRO via RS232 serial communication.

- Flat construction makes stepping on and off very easy.
- No need to turn it on and off - the scale does it for you.
- Low battery indicator for added convenience.

 **DETECTO**<sup>®</sup>  
A Division of Cardinal Scale Manufacturing Co.



## INTEGRAL CARRYING HANDLE

### 17-IN/32-CM EXTRA-WIDE WEIGHING PLATFORM FOR HEAVY PATIENTS

The SlimPRO's non-skid rubber mat allows for safe standing during measuring. The scale automatically turns on when a patient steps on the scale and turns off when they step off. Weight is conveniently held onscreen for 5 seconds after weighing.



### BATTERY POWER

The SlimPRO is powered by 4 AAA batteries, which are included with the scale. Auto-On and Auto-Off functions help conserve battery life.



### LARGE LCD READOUTS

The SlimPRO's big 1.5 in / 3.8 cm high LCD digits are easily viewable. Convenient units switching is performed through a switch on the back of the scale.



### SWIVEL FEET

DETECTO's SlimPRO can even weigh accurately on carpets. The special swivel feet allow the scale to be used on a variety of floor surfaces.




# LOW-PROFILE HEALTHCARE SCALES



Detecto  
Model SlimPRO  
Digital Low-Profile  
Floor Scale

<b>Capacity</b>	440 lb x 0.2 lb / 200 kg x 100 g
<b>Display</b>	Four digit, seven segment LCD display with 1.5 in / 3.8 cm high digits
<b>Display Area Size</b>	3.15 in W x 1.6 in H / 8 cm W x 4 cm H
<b>Overall Dimensions</b>	17 in W x 14.5 in D x 1.4 in H / 43 cm W x 37 cm D x 3 cm H
<b>Connectivity Port</b>	RS232 serial output (requires optionally purchased interface cable: Part No. 6600-1940)
<b>Power</b>	4 1.5V AAA alkaline batteries (included)
<b>Enclosure Construction</b>	Polycarbonate
<b>Units</b>	Kilograms, Pounds
<b>Features</b>	Auto-On / Auto-Off / Low Battery Indicator
<b>UPC Code</b>	809161189001
<b>Shipping Weight</b>	5 lb / 2.3 kg

*Detecto Scale reserves the right to improve, enhance, or modify features and specifications without prior notice.*

 **DETECTO**<sup>®</sup>  
A Division of Cardinal Scale Manufacturing Co.  
203 E. Daugherty, Webb City, MO 64870 USA  
Ph: 417-673-4631 or 1-800-641-2008 Fax: 417-673-5001  
[www.Detecto.com](http://www.Detecto.com)

**Sold By:**