



## TALKING SCALES (BILINGUAL ENGLISH/SPANISH)



## TALKING SCALES (ENGLISH AND SPANISH)

Easy to use for the vision-impaired

- Two capacities: 400 lb x 0.1 lb / 180 kg x 0.1 kg (SlimTalk)
   or 550 lb x 0.1 lb / 250 kg x 0.1 kg (SlimTalk XL)
- Extra-wide platform for patients of all sizes
- 5 carpet feet included
- Large 1.5 in/38 mm high LCD digits for easy reading
- 3 AAA batteries included
- Auto off function conserves battery life
- Clear, pleasant voice announcement of weight
  Very low-profile platform only 1.4 in/36 mm high for easy access

Sturdy construction for long-term

 Step-On/Off function for weighing without switching on scales


## NON-SLIP TEXU PLATFORM SURFA

durability



Slim TALK

Plastic platform with textured surface for weighing traction and user safety.



DETECTO.

Slim TALK

Plastic platform with black rubber mat for weighing traction and user safety.

	Model	SLIMTALK	SLIMTALKXL	
	Capacity	400 lb x 0.1 lb / 180 kg x 0.1 kg	550 lb x 0.1 lb / 250 kg x 0.1 kg	
	Languages		Two (English and Spanish)	
	<b>Overall Dimensions</b>	14.3 in W x 12.3 in D x 1.4 in H / 363 mm W x 312 mm D x 36 mm H		
ГCD	Display Size	3.2 in W x 1.7 in H / 81 mm W x 43 mm H		
Ц	Digit Height	1.5 in / 38 mm		
	Weight Display	Four digit, seven segment LCD with 1.5 in / 38 mm high digits		
Voice         English / Spanish / Off           Power         3 "AAA" batteries included (4.5 VDC)		English / S	panish / Off	
		included (4.5 VDC)		
	Indicators	lb, kg, Lo	lb, kg, Low Battery	
	Functions	Auto-ON (minimum weight of 11 lb / 5 kg is required)		
		Auto-OFF (10 seconds)		
Platform Construction         Plastic           Included         5 carpet feet, 3 AAA batteries, and product manual           Net Weight         5.7 lb / 2.6 kg		Plastic		
		eries, and product manual		
		/ 2.6 kg		
	Shipping Weight 7 lb / 3.2 kg		3.2 kg	
	UPC Code	809161302608	809161302707	



203 E. Daugherty, Webb City, MO 64870 USA Ph: 417-673-4631 or 1-800-641-2008 Fax: 417-673-2153

Mobile: m.detecto.com

DETECTO reserves the right to improve, enhance, or modify features and specifications without prior notice.

## SOLD BY:

© Copyright 2018 Cardinal Scale Mfg. Co. • Printed in USA • CAR/00/0218/C281