

STADIOMETERS

Mechanical, digital, and sonar stadiometers for precisely measuring height on all sizes of patients, ranging from infants up to adults



DETECTO

STADIOMETERS FOR ADULT AND PEDIATRIC PATIENTS

From the company renowned for weighing precision, DETECTO also manufactures a wide variety of stadiometers to select from based on your height measurement needs, ranging from wall-mount to freestanding, mechanical to digital and sonar. Baby length measuring stadiometers are available in both mechanical MLM and digital DLM types. Adult patients can be measured with the wall-mount mechanical 3PHTROD-WM or digital DHRWM, or the free-standing mechanical PHR or touchless sonar SONARIS.

MLM

Mechanical Baby Length Measuring



DETECTO's MLM stadiometer for infants measures from 2.5 in to 39 in / 65 mm to 1,000 mm in 0.125 in / 1 mm graduations. The durable, mechanical construction requires no electricity to operate. The calipers of the portable infant stadiometer conveniently fold in for easy transport or storage. It's light weight and easy to use. Soft ergonomic edges protect young patients from injury and maximize patient comfort while measuring.

DLM

Digital Baby Length Measuring



DETECTO's DLM stadiometer for infants measures from 13.75 to 31.5 in / 35 to 80 cm in 1/16 in / 1 mm graduations. Digital readouts are displayed on the easy-to-read LCD in feet/inches, inches, and centimeters. The DLM may be used stand-alone in the NICU or attached to the model MB130 pediatric scale. Powered by two 1.5V AAA batteries (included).

3PHTROD-WM

Adult Wall-Mount Mechanical Stadiometer



The 3PHTROD-WM measures from 30 to 78 in / 76 to 200 cm in 1/4 in / 1/2 cm graduations. Everything needed for operation is included with the stadiometer and installation takes seconds with the included hardware and built-in wall-mount bracket. Merely lift the spoon and request for the patient to step up to be measured. Easy to install and easy to operate.

PHR

Free-standing Mechanical Stadiometer

The PHR stadiometer measures from 4.5 in to 81 in / 11.5 cm to 205 cm in 0.125 in / 0.1 cm graduations visible on both sides of the height rod. There's no wall fastening necessary and it may be positioned against a wall or completely free standing. The PHR features a reversible head piece that lets you add a DETECTO DR400C floor scale at the base.



DHRWM

Adult Wall-Mount Digital Stadiometer

The DHRWM digital stadiometer measures from 43 in to 79 in / 110 cm to 200 cm in 0.5 in / 1 cm graduations. It's powered by 6 AA batteries (not included). Position the display where you need it with the wall-mount or desktop bracket. Crisp 1 in (25 mm) high LCD readouts and a 69 in (175 cm) cable length ensure ideal placement.



SONARIS

Free-standing or Wall-Mount Sonar Stadiometer

The SONARIS touchless sonar stadiometer measures from 0 to 9 ft / 0 to 275 cm in 0.2 in / 0.5 cm graduations. It provides instant digital readouts without pushing a button. It can be used stand-alone on a desktop/tabletop for mobility or wall mounted for permanent use. Powered by 6 AA batteries (not included).



Learn More:

<https://detecto.com/product/product-category/Stadiometers>

DETECTO reserves the right to improve, enhance, or modify features and specifications without prior notice.

SOLD BY:

DETECTO®

102 East Daugherty St., Webb City, MO 64870 USA
Toll-Free: (800) 641-2008 E-mail: detecto@cardet.com
www.Detecto.com