



- 192 oz x 0.05 oz / 5,500 g x 1 g
- Tare Button to Remove Dry Diaper Weight and Measure Fluid Output
- 80-hour Rechargeable Battery Pack with Automatic Shut Off Feature
- 5/8-in Deep Stainless Steel Removable Tray Aids Sample Retention and Allows Easy Clean Up
- 1.2-in/30-mm High Backlit LCD Easily Readable in Dimly-Lit Nurseries
- Scale is Completely Submersible Underwater





HYGIENIC AND HEAVY-DUTY

The watertight AC plug and stainless steel construction allow for easy clean up, prevent corrosion, and reduce the risk of bacteria harboring in nooks and crannies. The 5/8-in deep removable stainless steel tray with raised edges holds items neatly on the weighing platform and makes clean up effortless. Both of these features are a necessity for messy nursery applications.

DETECTO's Mariner* WPS12DT wet diaper scale meets all the hygiene demands of any active medical environment.

EASY-TO-USE, PORTABLE, AND CONVENIENT

The backlit blue LCD display offers large 1-in/25-mm high bold digits, readily visible in even the dimmest nursery conditions. The raised buttons allow for tactile functionality and easy operation. The UNIT key toggles through six different weight capacity display options in both decimals and fractions. The tare function conveniently allows the weight of a dry diaper to be removed from the total weight of a soiled diaper. The convenient 80-hour rechargeable battery pack with automatic shut-off feature conserves energy and allows the scale to be moved where it is needed the most.

Removable stainless steel tray makes sanitation effortless.	Backlit LCD is easily visible, even in darkened rooms so babies won't be disturbed.

Model	WPS12DT Wet Diaper Scales					
Capacity x Increments	12 lb x 0.1 oz 12 lb x 0.005 lb 12 lb x 1/8 oz	192 oz x 0.05 oz 192 oz x 1/8 oz 5,500 g x 1 g	Dimensions	Scale	7.5 in W x 9.8 in D x 2.3 in H / 190.5 mm W x 247.5 mm D x 57.0 mm H	
				Tray	9.0 in W x 6.0 in D x 1.1 in H / 228.6 mm W x 152.4 mm D x 27.9 mm H	
Zero	Established on power up routine and maintained by auto-zero circuitry		Tare		100% of full scale capacity	
Display / Digits	1.2-in / 30-mm hig (5) 1-in / 25-mm hi	•	Keypad Buttons		ON/OFF, UNITS/Backlight, SET, TARE	
Feet	6 rubberized feet for traction		Enclosure Rating		IP67 submersible	
Product Weight	4.9 lb / 2.2 kg		Operation Temperature		40° to 105°F (5° to 40°C)	
Shipping Weight	7.1 lb / 3.2 kg		Storage Temperature		40° to 140°F (5° to 60°C)	
Power	AC Adapter (with US plug included): 100-240 VAC input, 9V @ 1A output Battery (included): Ni-MH rechargeable 4.8 VDC 1800 mA (Note: a plug adapter is required for international use)		Please Note: The WPS12DT Battery Life		T Mariner® is not commercial dishwasher safe. 80 hours of continuous use with a fully recharged battery pack	
UPC Code	809161201901		Certification		NSF Approved	

Detecto Scale reserves the right to improve, enhance, or modify features and specifications without prior notice.

DETECTO. www.Detecto.com

102 E. Daugherty, Webb City, MO 64870 USA Ph: 417-673-4631 or 1-800-641-2008 Fax: 417-673-2153 www.Detecto.com

SOLD BY: