


solo[®] Physician Scale Parts Diagram



solo[®] Parts

Drawing Number	Part Number
1	7100-0218 Footing
2	7100-0219 Foot Switch
3	7100-0227 Left Cover for Load Cell (2 each)
4	7100-0228 Right Cover for Load Cell (2 each)
5	7100-0222 Column Base Cover
6	7100-0221 Friction Slider
7	7100-0224 Slider Assembly for Spoon
8	7100-0223 Spoon
9	7100-0225 Height Rod Top Cover
10	7100-0226 Collar for Height Rod
11	7100-0213 Cover for Serial Port
12	7100-0214 Cover for Batteries
13	7100-0215 Black Mat with Adhesive Backing
14	7100-0216 Platform Top Enclosure
15, 16, 17, 18, 19, 20, 21, 22, 2, 3	7100-0217 Scale Base Bottom Enclosure Assembly
23	7100-0211 Digital Indicator
24	7100-0210 Keypad Overlay
25	7100-0212 Cable to Indicator
26, 27, 28, 29	7100-0220 Wheel Kit
30	7100-0229 Indicator Mounting Bracket
AC ADAPTER	PD-AC

